
 <p>एफएसएसएआई fssai भारतीय खाद्य सुरक्षा और मानक प्राधिकरण Food Safety and Standards Authority of India स्वास्थ्य और परिवार कल्याण मंत्रालय Ministry of Health and Family Welfare</p>	Method for Determination of Iron in Fortified Rice Kernel by AAS		
Method No.	FSSAI.FRK.16.007.2023	Revision No. & Date	0.0
Scope	<p>The Scope of this Method is applicable for Quantification of Iron at 500 mg/Kg LOQ Level (with respect to the Sample) by using AAS.</p> <p>Limit of Detection 2.5 mg/Kg with respective to the Standard. Limit of Quantification 5.0 mg/Kg in with respective to the Standard. Limit of Quantification 500 mg/Kg in with respective to the Sample.</p>		
Safety & Precautions	<p>1. Concentrated Nitric Acid It is a Chemical which is corrosive to Metals. It causes severe skin burns and eye damage. It is toxic if inhaled. It is corrosive to the respiratory tract</p> <p>Following safety measures need to be taken during Handling of Concentrated Nitric Acid:</p> <ol style="list-style-type: none"> a) Do not breathe dust/fume/gas/mist/vapors/spray b) Wash face, hands and any exposed skin thoroughly after handling c) Wear protective gloves/protective clothing/eye protection/face protection d) Use only outdoors or in a well-ventilated area Keep away from heat/sparks/open flames/hot surfaces. e) No smoking f) Keep/Store away from clothing/ other combustible materials g) Take any precaution to avoid mixing with combustibles h) Keep only in original container i) Wear respiratory protection <p>2. Hydrogen Peroxide It is Oxidizing, Corrosive and Irritant chemical. Following safety measures need to be taken during Handling of Hydrogen Peroxide: When handling moderate-to-high concentrations of Hydrogen Peroxide in the workplace, ensure eyewash stations and safety showers are accessible, and use splash goggles, gloves, and an approved Vapor Respirator.</p>		
Principle	Weigh 0.50 g (\pm 0.05 g) of Grinded Sample Transfer to Microwave Digestion Cool Vessel. Add 2.0 mL Milli Q Water, 1.0 mL Hydrogen Peroxide, add 5 mL of Nitric Acid digest in microwave digester, extract the analyte in Nitric acid make up to 50 mL, Filter and Inject in AAS.		
Apparatus/Instruments	<ol style="list-style-type: none"> 1. Atomic Absorption Spectrometry (AAS) 2. Microwave Digester 3. Analytical Balance 4. Micro Pipettes (20 -200 μl) & (100 -1000 μl) <p>Note: The make & model of Instrument can be changed. However, the Instrument should be able to achieve the desired LOD value.</p>		
Materials and Reagents	<ol style="list-style-type: none"> 1. Concentrated Nitric Acid (Purity- 69%) 2. Hydrogen Peroxide (Purity -30%) 		

<p>Preparation of solutions</p>	<p>3. CRM Used: Iron</p> <p>A) PREPARATION OF INTERMEDIATE STOCK SOLUTION - 1 (100 mg/Kg)</p> <ol style="list-style-type: none"> 1. Transfer 10.0 ml from stock solution of iron (1000 mg/Kg) in 100 ml volumetric flask. 2. Add 5.0 ml nitric acid and made up the volume till 100 ml volumetric flask by Milli-Q water and mix by Vortex Shaker Mixer. <p>B) PREPARATION OF BRACKETING STANDARD SOLUTION (10 mg/Kg)</p> <ol style="list-style-type: none"> 1. Transfer 1.00 ml from Intermediate Standard Solution-1 of Iron (100 mg/Kg) in 10 ml volumetric flask. 2. Add 0.5 ml Nitric Acid and made up the volume till 10ml volumetric flask by Milli-Q water and mix by Vortex Shaker Mixer. <p>C) PREPARATION OF BLANK (5% NITRIC ACID)</p> <ol style="list-style-type: none"> 1. Transfer 7.25 mL of Nitric Acid (69%) in 100 mL Milli Q Water in Glass Bottle and Mix well. <p>D) PREPARATION OF CALIBRATION STANDARD SOLUTIONS</p> <ol style="list-style-type: none"> 1. Use Intermediate Standard Solution-1 for preparing Calibration Standard Solutions as mentioned in below Table. <table border="1" data-bbox="523 907 1503 1272"> <thead> <tr> <th>CAL. STANDARD SOLUTIONS</th> <th>ISS - 1 (mg/Kg)</th> <th>VOL. OF ISS - 1 (mL)</th> <th>VOL. OF NITRIC ACID (mL)</th> <th>FINAL VOL. (mL)</th> <th>FINAL CONC. (mg/Kg)</th> </tr> </thead> <tbody> <tr> <td>LS 6</td> <td>100</td> <td>8.00</td> <td>0.5</td> <td>10</td> <td>80.0</td> </tr> <tr> <td>LS 5</td> <td>100</td> <td>6.00</td> <td>0.5</td> <td>10</td> <td>60.0</td> </tr> <tr> <td>LS 4</td> <td>100</td> <td>4.00</td> <td>0.5</td> <td>10</td> <td>40.0</td> </tr> <tr> <td>LS 3</td> <td>100</td> <td>2.00</td> <td>0.5</td> <td>10</td> <td>20.0</td> </tr> <tr> <td>LS 2</td> <td>100</td> <td>1.00</td> <td>0.5</td> <td>10</td> <td>10.0</td> </tr> <tr> <td>LS 1</td> <td>100</td> <td>0.50</td> <td>0.5</td> <td>10</td> <td>5.0</td> </tr> </tbody> </table> <p>CAL : Calibration ISS : Intermediate Stock Solution VOL : Volume LS: Linearity Solution</p> <p>NOTE: Use freshly prepared Standard solutions for the analysis.</p>	CAL. STANDARD SOLUTIONS	ISS - 1 (mg/Kg)	VOL. OF ISS - 1 (mL)	VOL. OF NITRIC ACID (mL)	FINAL VOL. (mL)	FINAL CONC. (mg/Kg)	LS 6	100	8.00	0.5	10	80.0	LS 5	100	6.00	0.5	10	60.0	LS 4	100	4.00	0.5	10	40.0	LS 3	100	2.00	0.5	10	20.0	LS 2	100	1.00	0.5	10	10.0	LS 1	100	0.50	0.5	10	5.0
CAL. STANDARD SOLUTIONS	ISS - 1 (mg/Kg)	VOL. OF ISS - 1 (mL)	VOL. OF NITRIC ACID (mL)	FINAL VOL. (mL)	FINAL CONC. (mg/Kg)																																						
LS 6	100	8.00	0.5	10	80.0																																						
LS 5	100	6.00	0.5	10	60.0																																						
LS 4	100	4.00	0.5	10	40.0																																						
LS 3	100	2.00	0.5	10	20.0																																						
LS 2	100	1.00	0.5	10	10.0																																						
LS 1	100	0.50	0.5	10	5.0																																						
<p>Sample Preparation</p>	<p>PREPARATION OF SAMPLE SOLUTION</p> <ol style="list-style-type: none"> 1. Grind 50g sample as fine as possible. 2. Weigh 0.50 g (± 0.05 g) Grinded Sample. 3. Transfer to Microwave Digestion Closed (MDC) Vessel. 4. Heat Milli Q Water at 60 °C. 5. Add 2.0 mL of Hot Milli-Q water. 6. Add 1.0 mL Hydrogen Peroxide. 7. Add 5.0 mL of Nitric Acid. 8. Close the Microwave Vessel tightly. 9. Keep at Room Temperature for 5 minutes. 10. Keep the Vessel rotor in Microwave Digester. 11. Cool the Vessel at Room Temperature after Digestion. 																																										

	<p>12. Add 10 mL of Milli Q water.</p> <p>13. Mix well.</p> <p>14. Transfer to 50 mL Volumetric Flask.</p> <p>15. Volume make-up to 50 mL with Milli-Q water.</p> <p>16. Filter and use for the injection on AAS.</p>																										
Method of analysis	a) Instrument : AAS																										
	b) Equipment Conditions : As detailed in below Table																										
	<table border="1"> <tr> <td>Hallow Cathode Lamp</td> <td>Iron (as Fe)</td> </tr> <tr> <td>Lamp Current</td> <td>5 (mA)</td> </tr> <tr> <td>Absorption Wavelength</td> <td>372.0</td> </tr> <tr> <td>Slit Width(nm)</td> <td>0.2</td> </tr> <tr> <td>Signal - Type</td> <td>Atomic Absorption</td> </tr> <tr> <td>Signal Measurement</td> <td>Integration</td> </tr> <tr> <td>Oxidant</td> <td>Air</td> </tr> <tr> <td>Oxidant Flow(L/Min)</td> <td>13.5</td> </tr> <tr> <td>Acetylene Flow</td> <td>2</td> </tr> <tr> <td>Equation</td> <td>Linear</td> </tr> <tr> <td colspan="2" style="text-align: center;">Read Parameters</td> </tr> <tr> <td>Time(Sec)</td> <td>10</td> </tr> <tr> <td>Delay time(Sec)</td> <td>10</td> </tr> </table>	Hallow Cathode Lamp	Iron (as Fe)	Lamp Current	5 (mA)	Absorption Wavelength	372.0	Slit Width(nm)	0.2	Signal - Type	Atomic Absorption	Signal Measurement	Integration	Oxidant	Air	Oxidant Flow(L/Min)	13.5	Acetylene Flow	2	Equation	Linear	Read Parameters		Time(Sec)	10	Delay time(Sec)	10
	Hallow Cathode Lamp	Iron (as Fe)																									
	Lamp Current	5 (mA)																									
	Absorption Wavelength	372.0																									
	Slit Width(nm)	0.2																									
	Signal - Type	Atomic Absorption																									
	Signal Measurement	Integration																									
	Oxidant	Air																									
	Oxidant Flow(L/Min)	13.5																									
	Acetylene Flow	2																									
	Equation	Linear																									
	Read Parameters																										
	Time(Sec)	10																									
Delay time(Sec)	10																										
c) Microwave Digestion Program																											
<table border="1"> <thead> <tr> <th>SL. NO</th> <th>RAMPING STAGE</th> <th>HOLD TIME (Minutes)</th> <th>TEMP (°C)</th> <th>POWER (Watt)</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>01</td> <td>20</td> <td>180</td> <td>800</td> </tr> <tr> <td>2</td> <td>02</td> <td>10</td> <td>160</td> <td>800</td> </tr> <tr> <td>3</td> <td>03</td> <td>10</td> <td>140</td> <td>800</td> </tr> <tr> <td>4</td> <td>COOL DOWN</td> <td>10</td> <td>-</td> <td>-</td> </tr> </tbody> </table>	SL. NO	RAMPING STAGE	HOLD TIME (Minutes)	TEMP (°C)	POWER (Watt)	1	01	20	180	800	2	02	10	160	800	3	03	10	140	800	4	COOL DOWN	10	-	-		
SL. NO	RAMPING STAGE	HOLD TIME (Minutes)	TEMP (°C)	POWER (Watt)																							
1	01	20	180	800																							
2	02	10	160	800																							
3	03	10	140	800																							
4	COOL DOWN	10	-	-																							

Batch Organization	<p><u>Injection Sequence</u></p> <table border="1" data-bbox="580 248 1444 891"> <thead> <tr> <th>SL.NO.</th> <th>NAME OF INJECTIONS</th> <th>NUMBER OF INJECTIONS</th> </tr> </thead> <tbody> <tr><td>1</td><td>Blank</td><td>2</td></tr> <tr><td>2</td><td>Linearity Solution (LS) - 1</td><td>1</td></tr> <tr><td>3</td><td>Linearity Solution (LS) - 2</td><td>1</td></tr> <tr><td>4</td><td>Linearity Solution (LS) - 3</td><td>1</td></tr> <tr><td>5</td><td>Linearity Solution (LS) - 4</td><td>1</td></tr> <tr><td>6</td><td>Linearity Solution (LS) - 5</td><td>1</td></tr> <tr><td>7</td><td>Linearity Solution (LS) - 6</td><td>1</td></tr> <tr><td>9</td><td>Blank</td><td>2</td></tr> <tr><td>10</td><td>Sample Solution</td><td>1</td></tr> <tr><td>11</td><td>Blank</td><td>2</td></tr> <tr><td>12</td><td>Bracketing Standard Solution</td><td>1</td></tr> <tr> <td colspan="2">TOTAL INJECTIONS</td> <td>14</td> </tr> </tbody> </table>	SL.NO.	NAME OF INJECTIONS	NUMBER OF INJECTIONS	1	Blank	2	2	Linearity Solution (LS) - 1	1	3	Linearity Solution (LS) - 2	1	4	Linearity Solution (LS) - 3	1	5	Linearity Solution (LS) - 4	1	6	Linearity Solution (LS) - 5	1	7	Linearity Solution (LS) - 6	1	9	Blank	2	10	Sample Solution	1	11	Blank	2	12	Bracketing Standard Solution	1	TOTAL INJECTIONS		14
SL.NO.	NAME OF INJECTIONS	NUMBER OF INJECTIONS																																						
1	Blank	2																																						
2	Linearity Solution (LS) - 1	1																																						
3	Linearity Solution (LS) - 2	1																																						
4	Linearity Solution (LS) - 3	1																																						
5	Linearity Solution (LS) - 4	1																																						
6	Linearity Solution (LS) - 5	1																																						
7	Linearity Solution (LS) - 6	1																																						
9	Blank	2																																						
10	Sample Solution	1																																						
11	Blank	2																																						
12	Bracketing Standard Solution	1																																						
TOTAL INJECTIONS		14																																						
Calculation with units of expression	<p>a) Carry out a regression analysis and calculate Regression coefficient (R²) by analyzing the calibration standards by fitting the data into a linear regression curve, including zero as the response for the reagent blank.</p> $\text{Iron (mg/Kg)} = \frac{\text{Instrument Conc. (mg/Kg)} \times \text{Make-up Volume (mL)}}{\text{Sample Weight (g)}}$																																							
Results																																								
LOD & LOQ	<p>Limit of Detection 2.5 mg/Kg with respect to the Standard. Limit of Quantification 5.0 mg/Kg in with respect to the Standard. Limit of Quantification 500 mg/Kg in with respect to the Sample.</p>																																							
Reference	<p>RPT/MT/FRK/2023/001, Method Validation Report for Estimation of Iron in Fortified Rice Kernel by Using AAS.</p> <p>AOAC 2011.14: Determination of Minerals and Trace elements in Milk & Milk Products, Infant Formula, and Adult Nutrition.</p>																																							
Approved by	Scientific Panel on Methods of Sampling and Analysis																																							